

## How to add a case?

- 1) Go to "My cases" in the sidebar.
- 2) Click on "Add case" on the top right of the page.
- 3) Choose a name for your case and press "create"
- 4) Select the slides which you want to add to your case.
- 5) Click on "Add selected"

The slides now have been added to your case. To view your case, go back to the list by using the "back to list" button.

From:

<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:

[https://docs.pathomation.com/my/doku.php?id=how\\_to\\_add\\_cases&rev=1577356335](https://docs.pathomation.com/my/doku.php?id=how_to_add_cases&rev=1577356335)

Last update: **2019/12/26 13:32**

