


How to add a case?

1) Go to “My cases” in the sidebar.



2) Click on “Add case” on the top right of the page.

3) Choose a name for your case and press “Create”

4) Select the slides which you want to add to your case.

5) Click on “Add selected”

The slides now have been added to your case. To view your case, go back to the list by using the “back to list” button.

From:

<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:

https://docs.pathomation.com/my/doku.php?id=how_to_add_cases&rev=1577359268

Last update: **2019/12/26 14:21**

