



How to add cases?

1) Go to “My cases” in the sidebar.



2) Click on “+ New case” on the top right of the page.



3) Choose a name for your case and press “Create”

New case



* Name of new case

Cancel

Create

4) Select the slides which you want to add to your case.

5) Click on “>” button

Test Id: 744 ← Back to cases Preview

Slides Meta data

My Slides

- My Pathomation
 - bladder
 - lung
 - test
 - blur53.JPG
 - CMU-1.szi
 - CMU-3.jpg
 - Slide3.jp2

Test

<input type="checkbox"/>	Order	Thumbnail	Name	Path
 Select slides from tree and add them to case.				

The slides now have been added to your case. To view your cases, go back to the list by using the “Back to cases” button, or preview your newly created case by clicking on “Preview” button.

Test  Id: 744

[← Back to cases](#)


[Preview](#)

Slides Meta data

My Slides

- My Pathomation
 - bladder
 - lung
 - test
 - blur53.JPG
 - CMU-1.szi
 - CMU-3.jpg
 - Slide3.jp2**

[>](#)

<input type="checkbox"/>	Order	Thumbnail	Name	Path
<input type="checkbox"/>	▼ ▲		Slide3.jp2	My Pathomation/Slide3.jp2

From:
<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:
https://docs.pathomation.com/my/doku.php?id=how_to_add_cases&rev=1644842266

Last update: **2022/02/14 15:37**

