



## How to add cases?

A case can only be created from within a collection.

2) Click on “+ New case” on the top right of the page.



3) Choose a name for your case and press “Create”



4) Select the slides which you want to add to your case.

5) Click on “>” button


Test [↗](#) Id: 744 [← Back to cases](#) [Preview](#)

Slides    Meta data

My Slides

- My Pathomation
  - bladder
  - lung
  - test
    - blur53.JPG
    - CMU-1.szi
    - CMU-3.jpg
    - Slide3.jp2**

Test 🗑️

<input type="checkbox"/>	Order	Thumbnail	Name	Path
 Select slides from tree and add them to case.				

[➤](#)

The slides now have been added to your case. To view your cases, go back to the list by using the “Back to cases” button, or preview your newly created case by clicking on “Preview” button.


Test [↗](#) Id: 744 [← Back to cases](#) [Preview](#)

Slides    Meta data

My Slides

- My Pathomation
  - bladder
  - lung
  - test
    - blur53.JPG
    - CMU-1.szi
    - CMU-3.jpg
    - Slide3.jp2**

Test 🗑️

<input type="checkbox"/>	Order	Thumbnail	Name	Path
<input type="checkbox"/>	▼   ▲		Slide3.jp2	My Pathomation/Slide3.jp2

[➤](#)

From:  
<https://docs.pathomation.com/my/> - My Pathomation

Permanent link:  
[https://docs.pathomation.com/my/doku.php?id=how\\_to\\_add\\_cases&rev=1662383600](https://docs.pathomation.com/my/doku.php?id=how_to_add_cases&rev=1662383600)

Last update: **2022/09/05 16:13**

