



# How to Make a Course

The Courses option allows you to build content consisting of interactive images and text that can be shared with students or collaborators. You can begin creating a new course by clicking the “+New course” button (1) After the course is created you can add/edit the content by opening the edit panel (2) You can preview the course content you have created by clicking (3) You can share a course in the same way as you can images or collections by clicking (4).



Below is the editor window. Here you can add new content (1), which will prompt you to select which type of content you want to add (2), opting for Slide or Collection allows you to simply pick any of your available slides. You can configure what users following the course will be able to see via the viewport configuration button (3). Here you can choose what elements of the user interface (zoom, rotation etc.) and slide details (filename, barcode etc.) will displayed to students. You can preview the course with the preview button (4).



From:  
<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:  
[https://docs.pathomation.com/my/doku.php?id=how\\_to\\_make\\_a\\_course&rev=1699538907](https://docs.pathomation.com/my/doku.php?id=how_to_make_a_course&rev=1699538907)

Last update: **2023/11/09 17:08**

