



How to share a case?

1. Go to your cases via “My cases” on the sidebar.



2. All your cases are shown in the gallery. On every case, you can see the share button. Click share button and a new window pops up.



3. There are different ways to share your case.

1. You can copy the link and share it with others
2. You can embed your slide in your own web content
3. You can generate a QR-code for your slide
4. You can email the link of shared slide
5. You can directly send the link of shared slide via various social network service such as *Facebook, Twitter, WhatsApp, Messenger or LinkedIn.*

Also you can change the visibility option of your share. **Unlisted shares** can only be accessed by users with a valid link. **Public shares** can also be accessed by link, but are accessible from public listings too and anyone can view them.

Choose the options that suits your need best.

Share Test



Test

Created on Feb 10, 2022, 16:33

Contains 1 slide(s)

Share method

Visibility ⓘ: ☐ Public ☒ Unlisted

Link



Embed



QRcode



Email



Facebook



Twitter



WhatsApp



Messenger

ADVANCED (more options and properties)

Close

4. By clicking on “ADVANCED (more options and properties)” you can finetune your share. You can password-protect your share, set an expiration date for your share or modify viewport to your needs.

Advanced sharing

Test

Created on Feb 10, 2022, 16:33

Contains 1 slide(s)

Share method

Visibility ⓘ: ☐ Public ☒ Unlisted

Link

Embed

QRcode

Email

Facebook

Twitter

WhatsApp

Messenger

LinkedIn

Share options

Password-protect the share:

Password

Put an expiration on the share:

Select date

Viewport options

Channels

Default ☒

Default ☒

Time frame: 1 / 2

Z-stack: 2 / 2

Overview ⓘ: Visible ☒ Collapsed ☐

Dimensions control ⓘ: Visible ☒ Collapsed ☐

Barcode ⓘ: Visible ☒ Collapsed ☐

Rotation control ⓘ: Visible ☒ Collapsed ☐

Filename ⓘ: Visible ☐ Collapsed ☐

Scale line ⓘ: Visible ☒ Collapsed ☐

Snapshot button ⓘ: Visible ☐

Color adjustments ⓘ: Visible ☐

From:
<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:
https://docs.pathomation.com/my/doku.php?id=how_to_share_a_case&rev=1644841353

Last update: **2022/02/14 15:22**

