


How to share a slide?

1. Go to your slides via “My slides” on the sidebar.



2. All your slides are shown in the gallery. On every slide, you can see the share button. Click share button and a new window pops up.



3. There are different ways to share your slide.

1. You can copy the link and share it with others
2. You can embed your slide in your own web content
3. You can generate a QR-code for your slide
4. You can email the link of shared slide
5. You can directly send the link of shared slide via various social network service such as *Facebook, Twitter, WhatsApp, Messenger or LinkedIn.*

Also you can change the visibility option of your share. **Unlisted shares** can only be accessed by users with a valid link. **Public shares** can also be accessed by link, but are accessible from public listings too and anyone can view them.

Choose the options that suits your need best.

Share Slide3.jp2 X

 My Pathomation/Slide3.jp2

Share method Visibility: Public Unlisted


Link


Embed


QRcode


Email


Facebook


Twitter


WhatsApp


Messenger

[ADVANCED \(more options and properties\)](#) [Close](#)

4. By clicking on “ADVANCED (more options and properties)” you can finetune your share. You can password-protect your share, set an expiration date for your share or modify viewport to your needs.

Advanced sharing

 My Pathomation/Slide3.jp2

Share method Visibility: Public Unlisted


Link


Embed


QRcode


Email


Facebook


Twitter


WhatsApp


Messenger



LinkedIn

Share options

Password-protect the share:

Put an expiration on the share:

Viewport options



Options: Overview Visible Collapsed
 Dimensions control Visible Collapsed
 Barcode Visible Collapsed
 Rotation control Visible Collapsed
 Filename Visible Collapsed
 Scale line Visible Collapsed
 Snapshot button Visible
 Color adjustments Visible

From:
<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:
https://docs.pathomation.com/my/doku.php?id=how_to_share_a_slide&rev=1644842175

Last update: **2022/02/14 15:36**

