

Interaction Modes with Case Collections

There are different interaction modes which can be set per Case Collection during a training session. Users of a training session will only be able to see a Case Collection in the way designated by Project Owners/Supervisors.

The interaction modes [can be adjusted at any time before or during a session for all or just a subset of users](#).


This allows Supervisors/Project Owners to hide advanced content until Users have finished the preceding content, allow browsing of slide content without scoring, facilitate consensus scoring with or without Users being able to see each other's answers, and much more.

The possible interaction modes are listed below:

- Board - Collaboration mode is enabled with chat functionality
- Browse - only the slides are displayed
- Consensus view - the consensus view is available for everybody to see
- Hidden - the case collection is not visible
- Locked - the case collection is locked and scores cannot be edited
- Self review - the consensus score as well the user's own score is displayed
- Review - the consensus score as well as all user's individual scores are displayed
- Self test - a form is displayed to the participant for each case. The correct answers are immediately displayed upon submitting the form.
- Test active - a form is displayed to the participant for each case

Advanced interaction modes (for supervisors only)

- Clinical information edit: enter the clinical information data of this case
- Consensus score edit: the user is presented with a form where they can enter the consensus score values. All from data from all participants is shown.



The screenshot shows the PMA.control test environment interface. At the top, it says "PMA.control test environment" and "User" with a dropdown arrow and a help icon. Below that, the title "Training Session Test" is displayed. A message box says "Before starting, please take a few moments to complete a short questionnaire" with a "Test form Done ✓" button. Under "Case collections", there is a project description: "A project to show how easy it is to start a project!". Below this, there is a table of case collections:

Case Collection	Interaction Mode
Test	Locked
Test 3 0/9 cases evaluated	Take the test by scoring each case
Test 2	Browse the slides

From:
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:
https://docs.pathomation.com/pathotrainer/doku.php?id=interaction_modes&rev=1704287676

Last update: **2024/01/03 16:14**

