

Interaction Modes with Case Collections

There are different interaction modes which can be set per Case Collection during a training session. Users of a training session will only be able to see a Case Collection in the way designated by Project Owners/Supervisors.

The interaction modes [can be adjusted at any time before or during a session for all or just a subset of users](#).

This allows Supervisors/Project Owners to hide advanced content until Users have finished the preceding content, allow browsing of slide content without scoring, facilitate consensus scoring with or without Users being able to see each other's answers, and much more.



The possible interaction modes are listed below:

- Board - Collaboration mode is enabled with chat functionality
- Browse - only the slides are displayed
- Consensus view - the consensus view is available for everybody to see
- Hidden - the case collection is not visible
- Locked - the case collection is locked and scores cannot be edited
- Self review - the consensus score as well the user's own score is displayed
- Review - the consensus score as well as all user's individual scores are displayed
- Self test - a form is displayed to the participant for each case. The correct answers are immediately displayed upon submitting the form.
- Test active - a form is displayed to the participant for each case

Advanced interaction modes (for supervisors only)



- Clinical information edit: enter the clinical information data of this case
- Consensus score edit: the user is presented with a form where they can enter the consensus score values. All from data from all participants is shown.

For making a customized interaction mode, see following page [Custom Interaction Modes](#)

PMA.control test environment User  






Training Session Test

Before starting, please take a few moments to complete a short questionnaire

 Test form Done ✓ 

Case collections

A project to show how easy it is to start a project!

 Test Locked
 Test 3 0/9 cases evaluated Take the test by scoring each case 
 Test 2 Browse the slides 

From: <https://docs.pathomation.com/pathotrainer/> - Pathotrainer

Permanent link: https://docs.pathomation.com/pathotrainer/doku.php?id=interaction_modes&rev=1745565510

Last update: **2025/04/25 10:18**

