

Participant management

On this page a trainer/supervisor can select how the trainees can interact with the software. There are different [interaction modes](#)


[pma.control_interaction_modes.mp4](#)

This is performed through the Participant Management menu in the session overview menu. For details on where to find this see [the session page article](#).

Primarily you will interact with participants by changing the interaction mode of the slides/case collections and their associated forms.

In the screenshot below you can see how the page looks for the User (left) and the Project Owner (right) when various option

1. Locked - the content is not accessible
2. Self test - the number of cases the User has scored appears next to Test 3
3. Browse - the user can browse, but not score or see consensus/other user scores for the slides



From: <https://docs.pathomation.com/pathotrainer/> - Pathotrainer

Permanent link: https://docs.pathomation.com/pathotrainer/doku.php?id=participant_management&rev=1704279564

Last update: 2024/01/03 13:59

