

Participant management

On this page a trainer/supervisor can select how the trainees can interact with the software. There are different [interaction modes](#)

[pma.control_interaction_modes.mp4](#)

This is performed through the Participant Management menu in the session overview menu. For details on where to find this see [the session page article](#).

Primarily you will interact with participants by changing the interaction mode of the slides/case collections and their associated forms.

In the screenshot below you can see how the page looks for the User (left) and the Project Owner (right) when various interaction modes ([a full list of which can be found here](#)) have been selected

1. Locked - the content is not accessible
2. Self test - the number of cases the User has scored appears next to Test 3
3. Browse - the user can browse, but not score or see consensus/other user scores for the slides

After selecting the interaction mode(s) desired for each case collection, the Project Owner/Supervisor must select user(s) to apply this state to (4), and click Apply or Apply All (5)

The screenshot displays two side-by-side views of the PMA control test environment. The left view is for a 'User' and the right view is for a 'Project Owner'. The Project Owner view shows a 'Participants management' section with dropdown menus for selecting interaction modes (Locked, Self test, Browse) for different tests. Below this is a table of users with checkboxes for selection. Red arrows and numbers 1-5 indicate specific actions: 1 points to the 'Locked' mode for 'Test', 2 points to 'Self test' for 'Test 3', 3 points to 'Browse' for 'Test 2', 4 points to the 'User' checkbox in the table, and 5 points to the 'Apply all' button.

From:
<https://docs.pathomation.com/pathotrainer/> - Pathotrainer

Permanent link:
https://docs.pathomation.com/pathotrainer/doku.php?id=participant_management&rev=1704286955

Last update: 2024/01/03 16:02

