

# Participant management

On this page a trainer/supervisor can select how the trainees can interact with the software. There are different [interaction modes](#)

[pma.control\\_interaction\\_modes.mp4](#)


This is performed through the Participant Management menu in the session overview menu. For details on where to find this see [the session page article](#).

Primarily you will interact with participants by changing the interaction mode of the slides/case collections and their associated forms.

In the screenshot below you can see how the page looks for the User (left) and the Project Owner (right) when various interaction modes ( [a full list of which can be found here](#) ) have been selected

1. Locked - the content is not accessible
2. Self test - the number of cases the User has scored appears next to Test 3
3. Browse - the user can browse, but not score or see consensus/other user scores for the slides

After selecting the interaction mode(s) desired for each case collection, the Project Owner/Supervisor must select user(s) to apply this state to (4), and click Apply or Apply All (5)



From: <https://docs.pathomation.com/pathotrainer/> - Pathotrainer

Permanent link: [https://docs.pathomation.com/pathotrainer/doku.php?id=participant\\_management&rev=1704286955](https://docs.pathomation.com/pathotrainer/doku.php?id=participant_management&rev=1704286955)

Last update: 2024/01/03 16:02

