

Sharing Slide(s)


PIMSLS allows you to share content in a variety of ways. There is a dedicated group for sharing content on the Home tab of the ribbon:



From top left to bottom right the icons are:

1. Quick Share current folder
2. Quick Share current slide
3. Quick Share currently selected ROI / annotation
4. Quick Share current grid
5. Quick Share current ROIs/annotations visible in the grid
6. Quick Share barcode of current slide
7. Quick Share thumbnail of current slide
8. Custom share - this option shows you all the above possibilities and gives you extra customization options

Once you click of any of the Quick Share options you will be shown a popup window that allows you to copy the shared link to clipboard or to scan / download the QR-code



Custom Share

If you want more control over what and how you're sharing content, you can click on the final "Custom Share" button of the group. You could say that that's our "universal" share button.


It allows for further customization of your share link, including:

1. Setting a password to protect your shared content
2. Setting an expiry date (e.g. students can only access it for the duration of a test)
3. Include or exclude the barcode, annotations, and annotation labels from the shared link
4. Give the share a name
5. Add comments


The 'Share with others' dialog box is shown with the following options:

- Select what to share:** A list of sharing options with corresponding icons. To the right of each option is a preview image showing how the content will appear when shared. The options are: Folder (multiple slides), Slide (single slide), Partial slide (ROI), Grid of Slides, Grid of regions of interest, Barcode, and Thumbnail.
- Is password protected:** A checkbox that is currently unchecked.
- Expires:** A checkbox that is currently unchecked.
- Include:** Three checkboxes are checked: Barcode, Annotations, and Annotations Text.
- Name:** A text input field.
- Comment:** A text input field.
- Buttons:** 'Link preview' and 'Cancel' buttons are at the bottom right.


Password protected preview



Preview



And once you press the “Create Link” button it will generate the share link and show a popup window with option to “Copy to clipboard” or “Download QR-Code”



From:
<https://docs.pathomation.com/pimsls/3.0.0/> - **PIMS LS 3.0.0**

Permanent link:
https://docs.pathomation.com/pimsls/3.0.0/doku.php?id=slides_sharing&rev=1730905061

Last update: **2024/11/06 17:57**

