



AI Analysis


With the rapid growth of AI, PIMS LS has integrated some of the best AI algorithm to perform analysis on whole slide, selected area or current view. It gives pathologist detailed image analysis to diagnose things better.



Once user selected AI algorithm from the menu he is shown a pop box with details about algorithm, selected ROI, whole slide and also cost associated with performing the analysis. Its important in case of ROI that dimensions of it are with in selected AI algorithm boundaries else it would not allow to perform analysis. Once criteria is met "Run AI Analysis" button is enabled to perform analysis.



One the analysis is done, results are loaded in bottom section. It clearly indicates if analysis was successful or there was any error.



From the bottom panel user can click on info icon to get details on either success or failure. There is a side bar panel as well that shows up image results in case of success and we can navigate to different areas through it.

Viewport

AI

Default

Pre-set layout

CMU-1.mrx

100X

80X

40X

20X

10X

5X

2X

1X

1000 µm

AI analysis information for task 74f6fc12-d376-4852-8f4d-3b54659407f4

Task details

ID	74f6fc12-d376-4852-8f4d-3b54659407f4
Slide path	_sys_ref/annotations/CMU-1.mrx
Status	Failed
Request submitted on	5/12/2025, 4:42:07 PM
Request completed on	5/12/2025, 4:42:14 PM

AI algorithm details

ID	mindpeak-roi-k67
Name	AI Tumor Detection Service ROI
Description	Detect and quantify KI-67 within predefined RoIs of breast cancer tissue
Type	TumorDetection
Method	KI67
Specimen source site	Breast
Vendor	Mindpeak
Regulatory approval	EU: CE-IVD Outside the EU: Research Use Only

Error

Error codes	N/A
-------------	-----

Close

Completed out of 1 Total

CMU-1.mrx

AI

From:

<https://docs.pathomation.com/pimsls/3.1.0/> - **PIMS LS 3.1.0**

Permanent link:

<https://docs.pathomation.com/pimsls/3.1.0/doku.php?id=aianalysis>

Last update:

2025/05/12 20:01

