

# Working with grids

PMA.studio uses viewport as a grid system. UI provides a separate tab “GRID” to perform grid operations on the viewport.



“GRID” tab consists of 6 sections:

1. Size
2. Synchronize
3. R/C
4. Slide
5. Snapshot
6. Share


## Size

Size section allows you to break viewport as a grid system with different sizes. This includes grid sizes of 1x1,2x2,3x3 and even you can customized the grid size according to your need.




### (1x1 Grid size)

First grid option allows you to divide the viewport into 1x1 column grid.



### (2x2 Grid size)

Second grid option divides the viewport into 2x2 column grid.



### (3x3 Grid size)

Third grid option divides the viewport into 3x3 column grid.



### (Custom Grid size)

You can even customize the grid according to your need by using the custom grid option.



### Synchronize

This section allows you to perform different actions as a whole on individual grid column or as a whole grid system. It is further divided into 4 sub actions:



- 1. Toggle synchronous navigation between cells in the current grid

This button allows you to perform set of synchronous actions on whole grid.

- 1. Clear all grid cells


This button allows clearing all the grid cells at the same time

- 1. Clear current cell only

This button allows to clear the selected grid cell only

- 1. Open current slide into separate viewport panel

This button opens the current selected slide into a separate viewport tab



From:  
<https://docs.pathomation.com/pimsls/3.1.1/> - PIMS LS 3.1.1

Permanent link:  
[https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=slides\\_grid&rev=1650996405](https://docs.pathomation.com/pimsls/3.1.1/doku.php?id=slides_grid&rev=1650996405)

Last update: 2022/04/26 21:06

