

Meta-data

With meta-data, we identify all data that surrounds a slide, but that is not graphical. The actual pixels (or viewport representation) of a slide are *not* meta-data. [Annotations](#) can be considered meta-data in the semantic sense, but they're link to graphical activities as well. So [annotations are discussed elsewhere in this manual](#).

There is [a blog article](#) that serves as a complementary tutorial to this blog entry. In [the article](#), we consider 3 types of meta-data in PMA.core:

Meta-data	Location	Storage
Intrinsic meta-data	Embedded in the slide's file format	File system or cloud FS
User-captured meta-data	Stored in PMA.core	PMA.core RDBMS (SQLServer)
External meta-data	Stored in own (customer) database	SQLServer or Oracle

Intrinsic meta-data

This type of data can be consulted via the Image Management's [left panel](#), or more directly via the [API](#).

User-captured meta-data

In the [forms management section](#), data input forms can be designed that are subsequently assigned to virtual slides. These forms can then be filled in by the system's users while viewing slides. The data that users submit for each form are also available within this section.

External meta-data

From:
<https://docs.pathomation.com/pma.core/2.0.1/> - **PMA.core 2.x**

Permanent link:
https://docs.pathomation.com/pma.core/2.0.1/doku.php?id=meta_data&rev=1648460179

Last update: **2022/03/28 12:36**

