


Heatmap statistics

Heatmap of a slide is a calculated statistic by PMA.core that shows the regions that users request more often in a slide and thus be of more importance. This is displayed in a colorful image with more reddish colors indicating regions of the image requested more times by users, as is the case with typical heatmaps in other fields. At the bottom of the heatmap you can see the color legend describing what each color indicates the number of times a region was requested by users

To generate a Heatmap for a specific slide you need to:

1. Select a date range to request usage data. You can narrow this range to your needs
2. Select the folder that contains the slide you want in the middle column
3. Select the slide you want in the last column
4. After the thumbnail of the slide is shown click on "Calculate Heatmap"



From:

<https://docs.pathomation.com/pma.core/2.0.2/> - **PMA.core 2.x**

Permanent link:

https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=report_heatmap&rev=1648648127

Last update: **2022/03/30 16:48**

