


Security

Security is increasingly important. As PMA.core has been deployed in increasingly complex scenarios over the years, its security features have evolved, too.

Security pertaining to root-directories is situated at two levels:

- Security features that enable root-directories to access content, such as:
 - Configure public/secret key combinations for S3 resources
 - Configure account credentials to be used when accessing a UNC network resource path
- Prevent users from access mounted content through root directories that they are or are not allowed to do
 - Define Access control lists



Add mounting point

The network name cannot be found

Instance: Default

Path: c:\pathomation\secure_ws

Impersonation username: pma_read

Impersonation password: [masked]


So when defining local hard disk paths, make sure the impersonation options are left blank.

Network storage (UNC paths)


Pathomation runs under a certain application pool. This application pool is associated with a user identify, which may not have access to the network path that you try to access. Giving access for the application pool to access the network resource may be difficult for a variety of reasons.

If you can't immediately access the network path with default (i.e. application pool) credentials, you can provide additional information.

In the case below we've created a dedicated pma_read user that is permitted to access the shared [\\MALTA1767\reference](#) path:



We can enter this as a path for the mounting point, and add the impersonation information for our pma_read user:




Path: \\MALTA1767\reference

Impersonation username: pma_read

Impersonation password: [masked]

Impersonation domain name: MALTA1767

The mounting point shows up, and you can activate the View slides tab to inspect its content:



If the credentials are faulty, an error appears

S3 storage

Azure storage

Public vs private

As you have more users and more root-directories, it becomes undesirable that everybody is allowed to see everything.


Therefore, root-directories can be marked “public” or “private”:

Public root directories are marked “public”, it means every user has access to them. They can be accessed by anybody who is a registered user in the PMA.core user repository.

Private root directories are marked “private”, it means only select users can see the content. They are only accessible by those who have been explicitly given access to be allowed to access the folder through the directory's [access control list](#).

Access control lists

Once marked private, you can select what users are allowed to see the content of the root directory, and which ones aren't: Do this by pressing the “Edit access control list” link after you selected the “private” option:



Last update: 2022/02/11
13:09

rootdir_security https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=rootdir_security&rev=1644574185

From:

<https://docs.pathomation.com/pma.core/2.0.2/> - **PMA.core 2.x**

Permanent link:

https://docs.pathomation.com/pma.core/2.0.2/doku.php?id=rootdir_security&rev=1644574185

Last update: **2022/02/11 13:09**

