

Working with grids

PMA.studio uses viewport as a grid system. UI provides a separate tab “GRID” to perform grid operations on the viewport.




“GRID” tab consists of 6 sections:

1. Size
2. Synchronize
3. R/C
4. Slide
5. Snapshot
6. Share


Size

Size section allows you to break viewport as a grid system with different sizes. This includes grid sizes of 1x1,2x2,3x3 and even you can customized the grid size according to your need.



(1x1 Grid size)

First grid option allows you to divide the viewport into 1x1 column grid.



(2x2 Grid size)

Second grid option divides the viewport into 2x2 column grid.



(3x3 Grid size)

Third grid option divides the viewport into 3x3 column grid.



(Custom Grid size)

You can even customize the grid according to your need by using the custom grid option.



Synchronize

From:

<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:

https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=slides_grid&rev=1650995506

Last update: **2022/04/26 20:51**

