

# PMA.studio

Welcome to PMA.studio, an integrated web-based environment for whole slide image (WSI) data, digital pathology, and virtual microscopy alike, developed by [Pathomation](#).

With PMA.studio pathomation has developed a true cockpit for imaging scientists and histopathologists.

the user is not locked into a fixed user interface, customization is a key characteristic of PMA-Studio, and a philosophy underpinning product development.


Guided by user experience and feedback, PMA.studio provides standard features which are considered essential and that should not be changed, but apart from those, the [user interface](#) can be extensively customized.

Whether you're looking to annotate slides, integrate distributed heterogeneous databases with your slides, or are looking to host collaborative long-distance tele-conferences, Pathomation's PMA.studio is a true cockpit for the pathologist and (imaging) research scientist alike.

PMA.studio features is built on top of PMA.view, and hosts a slew of innovative features:


## Annotations

PMA.studio comes with an Annotations panel where you can review and manage annotations for specific slides, You can draw your own annotations of different shapes. You can then export your annotations to a number of popula formats.



## Data management

Digital pathology involves more than just whole slide images. Oftentimes supplemental data is available in external data stores. These data stores can be consulted in addition to the image data itself. We call it “integrated pathology”



## Virtual slide tray


PMA.Studio offers to organize slides in a tray - it's a virtual shell where you can put specific slides or folders. Virtual tray allows you to review slide info - dimensions, file size, # of layers and channels, barcode, filled forms, as well as external information.

A screenshot of the PMA.Studio virtual slide tray. It shows a table with the following columns: Name, File name, Width, Height, Max Magnification, and p Per Pixel. The table contains several rows of slide information, including file names like "\_job\_information\the path\programs\FacCK...".


Name	File name	Width	Height	Max Magnification	p Per Pixel
PMA Case	_job_information\the path\programs\FacCK... C06_0000\F020000.mrx	8312	10624	40x	0.2427073
PMA Case	cases_autoload\06420467B_PDL_1_06420467B0.mrx	8094	10192	40x	0.2430019
PMA Case	cases_autoload\06420467B_PDL_1_06420467B0.mrx	8076	10208	40x	0.2427073
PMA Case	cases_autoload\06420467B_PDL_1_06420467B0.mrx	8476	10212	40x	0.2420019
PMA Case	cases_autoload\06420467B_PDL_1_06420467B0.mrx	8028	10438	40x	0.2430019
PMA Case	cases_autoload\06420467B_PDL_1_06420467B0.mrx	8094	10204	40x	0.2430019
PMA Case	cases_autoload\06420467B_PDL_1_06420467B0.mrx	8724	10208	40x	0.2420019
PMA Case	cases_autoload\06420467B_PDL_1_06420467B0.mrx	8076	17502	40x	0.2427073
PMA Case	cases_autoload\06420467B_PDL_1_06420467B0.mrx	8028	10438	40x	0.2427073
PMA Case					
PMA Case					
PMA Case					
PMA Case					
PMA Case					
PMA Case					

## Conferencing

Host or participate in virtual conferences that include joint and simultaneous exploring slides or grid of slides and annotating slides. Conferences can be secured with a password and allow for bi-directional synchronization.



## PMA.Studio Environment



From:  
<https://docs.pathomation.com/pma.studio/2.1.0/> - **PMA.studio 2.1.0**

Permanent link:  
<https://docs.pathomation.com/pma.studio/2.1.0/doku.php?id=start&rev=1680714777>

Last update: **2023/04/05 20:12**

