


The site manager aims to make it easy for end users to store connection data related to the different Pathomation server instances (PMA.core, Pathomation cloud...).

It can be initiated via :


- click on button “Site Manager” (first to the left on the icons toolbar)



- click on Menu item “Site Manager” on Menu “File”.



- Keyboard shortcut “Ctrl + M”.



It enables end users to organize connection data into folders, sites, bookmarks in addition to a couple of others functionalities (rename, duplicate...) The data generated by the site manager is stored inside file “PMA.Transfer\_data.json” on user’s default Windows directory :

| Name                     | Date modified        | Type        | Size |
|--------------------------|----------------------|-------------|------|
| ★ Favorites              | 12/14/2018 11:42 ... | File folder |      |
| 🔗 Links                  | 12/14/2018 11:42 ... | File folder |      |
| 🎵 Music                  | 12/14/2018 11:42 ... | File folder |      |
| ☁ OneDrive               | 1/2/2019 1:40 AM     | File folder |      |
| 🖼 Pictures               | 12/14/2018 11:42 ... | File folder |      |
| 🎮 Saved Games            | 12/14/2018 11:42 ... | File folder |      |
| 🔍 Searches               | 12/14/2018 11:42 ... | File folder |      |
| 📄 TwixTel                | 3/22/2018 12:45 PM   | File folder |      |
| 📺 Videos                 | 12/14/2018 11:42 ... | File folder |      |
| 📄 PMA.Transfer_data.json | 1/7/2019 1:30 AM     | JSON File   | 1 KB |

The content of this file should not be modified!

From:  
<https://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link:  
[https://docs.pathomation.com/pma.transfer/doku.php?id=about\\_site\\_manager&rev=1589595288](https://docs.pathomation.com/pma.transfer/doku.php?id=about_site_manager&rev=1589595288)

Last update: **2020/05/16 05:14**

