

What's new in version 3.0.0?

Pathotrainer 2.2.0 is a release of the PathoTrainer codebase.

It is a minor release that contains several new features helpful for multi-stage training programs, extracting actionable information from global evaluation programs, and reducing the admin burden for these.

New Features

- If you have more than one PMA.core server with us e.g. in USA, LATAM, EU, Asia, Africa and Oceania, it is now possible to turn on automatic server assignment in a training session so that users will be connected to the fastest available server.
 - *It is still possible to manually assign a dedicated server per training session*
- Country tagging for users - Users can now input or be assigned to a country
- Training Session results are now searchable/filterable by country
- Users can be searched for by country when adding participants to a training session
- Progression based unlocking - It is now possible to make Case Collections in a Training Session unlock sequentially as the Trainee completes them

Older Versions

What's new in version 2.3.1?

Pathotrainer 2.3.1 is a release of the PathoTrainer codebase.

It is a patch release that resolved behaviour errors and performance issues with multiperson trainings that include large data sets AI generated annotations.

Bug fixes

- PIMS-4164: Temporary annotations are no longer saved but deleted at the end of a session
- PIMS-4044: Case attachments moved to the recycle bin become unavailable for users
- PIMS-4026: Erased Annotations are permanently deleted and cannot reappear
- PIMS-4023: Infinite Loading State When Switching Cases in Live Training Session resolved
- PIMS-4018: Tumor Complete Strong annotation count updates after adding new annotation
- PIMS-2609: Annotation toolbar is enabled in course mode
- PIMS-2484: Add option to control if scoring a case should advance the user to the next case automatically
- PIMS-1135: Keyboard shortcuts explanations language errors in Dutch corrected
- PIMS-1133: Saving consensus scores timeout error fixed
- PIMS-1131: Importing cases: Checked cases that are not on the current page of results when you click 'Import' are not imported.
- PIMS-910: when using a laptop, the UI scales properly and the user is able to access all the case

navigation and Live Mode buttons

What's new in version 2.3.0?

Pathotrainer 2.3.0 is a release of the PathoTrainer codebase.

It is a minor release that added a new feature to include AI generated annotations and results in images and incorporate these in the learning path of the trainees.

New Features

- Show prepopulated AI annotations + calculated score results in the PathoTrainer UI to enable trainees and trainers to leverage AI-assisted diagnostic insights during training sessions, enhancing learning outcomes and providing real-world AI workflow experience

Bug fixes

- PIMS-3912: Annotations shift position (offset to left) during zoom actions in Slide View.
- PIMS-3904: AI logo/panel slow performance results in wrong slides/case display mismatch
- PIMS-3903: Participants of training sessions with 'trainee' role cannot see the Mindpeak annotations changes in live sessions

What's new in version 2.2.0?

Pathotrainer 2.2.0 is a release of the PathoTrainer codebase.

It is a minor release that contains several new features helpful for multi-stage training programs, extracting actionable information from global evaluation programs, and reducing the admin burden for these.

New Features

- If you have more than one PMA.core server with us e.g. in USA, LATAM, EU, Asia, Africa and Oceania, it is now possible to turn on automatic server assignment in a training session so that users will be connected to the fastest available server.
 - *It is still possible to manually assign a dedicated server per training session*
- Country tagging for users - Users can now input or be assigned to a country
- Training Session results are now searchable/filterable by country
- Users can be searched for by country when adding participants to a training session
- Progression based unlocking - It is now possible to make Case Collections in a Training Session unlock sequentially as the Trainee completes them

Bug fixes

- PIMS-2751: Overlapping due to long text fields was fixed: User dropdown (logout etc.), Case name between the previous and next buttons, upper left side of page session name does not cover part of the WSI, Filename is displayed as seen in the image on smaller screens in order to avoid overlapping on the fullscreen button.

What's new in version 2.1.0?

Pathotrainer 2.1.0 is a release of the PathoTrainer codebase.

It is a minor release that expands on the annotation capabilities to allow for 'teacher' annotations.

New Features

- When using the consensus user edit mode, you are able to create 'Supervisor' level annotations that only reveal to the Trainee upon submission of their own answer.

What's new in version 2.0.0?

Pathotrainer 2.0.0 is a release of the Pathotrainer codebase (previously known as PMA.control).

It is a major release that contains multiple new features including:

- Trainees can now create annotations visible only to Supervisors within a Training Session,
- Admins can now create always-online self-study courses complete with automatic certification, quizzes,
- and much more - see the full notes below:

New Features

- CTRL-27: Made it possible to clone / duplicate forms
- CTRL-48: Enabled measurements in the training session slide viewer
- CTRL-60: Added new data export options including csv, normalized data
- CTRL-63: Documentation tab now directs users to the Pathotrainer wiki
- CTRL-64: A certification tab has been added to allow automatic certification for training sessions
- CTRL-76: API calls can now optionally filter by session name
- CTRL-82: Added shortcuts to make it easier to manage a project
- CTRL-85: Integrated changes to ensure compatibility with PMA.core 3
- CTRL-91: A pop up message has been made to increase awareness of the "Sync" function.
- CTRL-93: Course Mode - functionality to support the creation of standalone self-learning resources has been added.
- CTRL-94: The functionality to make Quizzes has been added.
- CTRL-95: Additional zoom level controls have been added at the bottom of the viewer UI
- CTRL-96: Users can now make annotations in a session outside of Live Mode. These are only visible to supervisors, and the supervisor is able to filter annotations by user(s)
- CTRL-97: Sharing of WSIs and ROIs is now possible

- CTRL-98: Project owners / Supervisors can now set how all image labels will be displayed (i.e. filenames | captions | no label) with a single click in the Training Session options.

Bug fixes

- CTRL-13: The 'Reset Score' button no longer appears when resetting the score is not permitted.
- CTRL-16: Fixed an issue that would allow scoring on the wrong case
- CTRL-17: Fixed an issue that could transfer form data from one case to another in rare circumstances
- CTRL-23: Fixed an issue causing some scores to not appear correctly in training sessions
- CTRL-51: Fixed the way the "Insert link" button was displayed in training sessions
- CTRL-88: Changed License expiry notification visibility to only admin users
- CTRL-145: Fixed bug in which WSI content would not load in Course Mode

From:

<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:

https://docs.pathomation.com/pathotrainer/doku.php?id=release_notes&rev=1788784265

Last update: **2026/09/07 14:31**

