

# Context menu for folders

## Local site

|   | Slidename           | Slidesize | FileType    | Last Modified |
|---|---------------------|-----------|-------------|---------------|
| ↑ | ..                  |           |             |               |
| 📁 | Intel               |           | File folder |               |
| 📁 | local samples       |           | File folder |               |
| 📁 | PerfLogs            |           | File folder |               |
| 📁 | Program Files       |           | File folder |               |
| 📁 | Program Files (x86) |           | File folder |               |
| 📁 | scripts             |           | File folder |               |
| 📁 | Shared folder       |           | File folder |               |
| 📁 | Slide Samples       |           | File folder |               |
| 📁 | SvmCache            |           | File folder |               |

Upload folder to remote site

Initiated via mouse right click within the table view borders, it's composed of one menu item : "Upload folder to remote site".

Upload folder to remote site :

Has been be discussed in details earlier on this tutorial.

## Remote site

No context menu options available yet for folders on remote site connections.

From: <http://docs.pathomation.com/pma.transfer/> - **PMA.transfer**

Permanent link: [http://docs.pathomation.com/pma.transfer/doku.php?id=context\\_menu\\_for\\_folders](http://docs.pathomation.com/pma.transfer/doku.php?id=context_menu_for_folders)

Last update: **2021/07/08 15:51**

