



Windows

- Extract file PMA.transfer.zip into a destination a folder
- Navigate to destination folder

| Name | Date modified | Type | Size |
|---|--------------------|------------------------|----------|
|  PMA.transfer-2.0.3.72.exe | 1/24/2020 11:08 AM | Application | 9,483 KB |
|  PMA.transfer.ini | 3/17/2019 10:11 PM | Configuration settings | 1 KB |

- Double click on file PMA.transfer-2.0.3.x.exe

The Zip file is optionally shipped with a PMA.transfer.ini file containing a preconfigured connection data, it's loaded automatically when the app is launched.

MacOS

From:
<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:
<https://docs.pathomation.com/PMA.transfer/doku.php?id=installation&rev=1598914769>

Last update: **2020/09/01 01:59**

