


The node "My Sites" is the root of the site manager tree, it cannot be deleted or renamed.



- Click on button "New Site" creates a new site inside "My Sites".
- Click on button "New Folder" creates a new folder inside "My Sites".

From:  
<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:  
[https://docs.pathomation.com/PMA.transfer/doku.php?id=my\\_sites&rev=1589595289](https://docs.pathomation.com/PMA.transfer/doku.php?id=my_sites&rev=1589595289)

Last update: **2020/05/16 05:14**

