

# Release history

| Version                            | Revision   | Released       |
|------------------------------------|------------|----------------|
| PMA.transfer 2.0.0                 | 50         | May 2019       |
| <a href="#">PMA.transfer 2.0.1</a> | 55         | July 2019      |
| <a href="#">PMA.transfer 2.0.2</a> | 60         | September 2019 |
| <a href="#">PMA.transfer 2.0.3</a> | 72         | January 2020   |
| <a href="#">PMA.transfer 2.1.0</a> | 85         | May 2020       |
| <a href="#">PMA.transfer 2.1.1</a> | 87         | June 2020      |
| <a href="#">PMA.transfer 2.2.0</a> | 103        | August 2020    |
| <a href="#">PMA.transfer 2.2.1</a> | 111        | February 2021  |
| <a href="#">PMA.transfer 2.2.2</a> | 114        | July 2021      |
| <a href="#">PMA.transfer 2.2.3</a> | 9d6e296f92 | August 2021    |
| <a href="#">PMA.transfer 2.2.4</a> | d67f7abd23 | February 2021  |

## What's new in version 2.0.1

### General improvements :

- Log file menu item to open directly the current log file in addition to log files' location
- Added file format to “inspect slide” interface

### Bug fixes :

- VSI slides' transfer issue
- Dealing with leading/trailing spaces in host & username text fields and trailing slashes for host as well

## What's new in version 2.0.2

### General improvements :

- Ask for confirmation before closing the app
- Auto detect white or black text color depending on text area background color
- Add recursive folder download feature

### Bug fixes :

- Fix thumbnail URL encoding issue on local site
- Tests against all supported file formats, bug fixing

## What's new in version 2.0.3

### New features :

- Bypass proxy
- Delete directories on server side (only if empty and not a root one)

- Rename directories on server side (only if empty and not a root one)
- Delete slides on server side
- Rename slides on server side
- Add copy full path option both locally and remotely

#### General improvements :

- Directly download to selected folder on local site
- Keyboard input interaction
- New icon for root directories
- Better transfer speed metrics
- Better error messages
- Enhanced integrity check
- Adding counters to all bottom bar buttons

#### Bug fixes :

- Refresh table view for local site following successful downloads

#### **What's new in version 2.1.0**

#### New features :

- Interfacing with the new My Pathomation API
- Adding create/rename/delete folder features for My Pathomation accounts
- Adding rename/delete slide features for My Pathomation accounts
- Adding drag & drop functionality for download
- Release of a self-containing package for Windows (No installation of JRE required)

#### General improvements :

- Detection of PMA.core version on connection
- Harmonization of upload/download UI/UX
- Cleaning up UI for slides without barcode
- Defining explicit location while adding new folders

#### Bug fixes :

- Disabling of upload/download buttons on disconnection

#### **What's new in version 2.1.1**

#### General improvements :

- Compliance with PMA.start and PMA.start dotnetcore on Windows
- Adding support of VPN connections

#### **What's new in version 2.2.0**

### New features :

- Server to server transfers
- Multiple selection of slides for transfer

### General improvements :

- Compliance with PMA.start dotnetcore on MacOS & Unix
- Improved error messages (transfer of empty folders, transfers to location that require elevation of privileges..)
- Improved integrity check process

### Bug fixes :

- Better handling of failed connections
- Fixes on the site manager

## **What's new in version 2.2.1**

### New features :

- User quota enforcement

### General improvements :

- Error handling for uploads to Azure directories
- Improved error handling for server to server transfers
- Support for My Pathomation "add slide" API endpoint & "caller ID" parameter for slides uploaded via PMA.transfer
- Support for Keyboard key <Enter>

### Bug fixes :

- Fixes on the server to server progress bar


## **What's new in version 2.2.2**

### New features :

- Improved reporting of successful slide transfers (interactive vs. console)

### Bug fixes :

- Slide (file)size is reported correctly again
- Connectivity to My Pathomation is operational again



From:

<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/PMA.transfer/doku.php?id=release\\_history&rev=1677499209](https://docs.pathomation.com/PMA.transfer/doku.php?id=release_history&rev=1677499209)

Last update: **2023/02/27 15:00**

