


# PMA.transfer wiki start page

PMA.transfer is a slide transfer desktop application aimed to make slide transfer an easier process by remedying to a couple of frequent recurring problems (transfer failure, timeouts, Firewall issues, security, integrity check...)



Pathomation makes the PMA.transfer program available free of charge.

The software offers the following benefits:

- Knowledge of what a slide actually is: you select a .mrxs file; the application automatically figures out to include the other files in the process.
- Runs over HTTP(S): no need to punch through firewalls.
- Error recovery in case the network connection is interrupted.
- Integrity check to verify the transferred files are not corrupt.

From: <https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link: <https://docs.pathomation.com/PMA.transfer/doku.php?id=start&rev=1624294785>

Last update: **2021/06/21 19:59**

