

PMA.transfer is available for Windows, MacOS and Unix. The minimal distribution required for each is :

- Windows 7 or later
- MacOS 10.14 or later

It's worth mentioning that the Java runtime environment is no longer required as it is included out of the box in the self-containing package for each platform : Windows, MacOS and Unix.

From:

<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/PMA.transfer/doku.php?id=system\\_requirements&rev=1598914730](https://docs.pathomation.com/PMA.transfer/doku.php?id=system_requirements&rev=1598914730)

Last update: **2020/09/01 01:58**

