

### General improvements :

- Log file menu item to open directly the current log file in addition to log files' location
- Added file format to “inspect slide” interface

### Bug fixes :

- VSI slides' transfer issue
- Dealing with leading/trailing spaces in host & username text fields and trailing slashes for host as well

From:

<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.0.1&rev=1589519799](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.0.1&rev=1589519799)

Last update: **2020/05/16 05:14**

