

### General improvements :

- Ask for confirmation before closing the app
- Auto detect white or black text color depending on text area background color
- Add recursive folder download feature

### Bug fixes :

- Fix thumbnail URL encoding issue on local site
- Tests against all supported file formats, bug fixing

From:

<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.0.2&rev=1589595290](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.0.2&rev=1589595290)

Last update: **2020/05/16 05:14**

