

# What's new in version 2.1.0

## New features :

- Interfacing with the new My Pathomation API
- Adding create/rename/delete folder features for My Pathomation accounts
- Adding rename/delete slide features for My Pathomation accounts
- Adding drag & drop functionality for download
- Release of a self-containing package for Windows (No installation of JRE required)

## General improvements :

- Detection of PMA.core version on connection
- Harmonization of upload/download UI/UX
- Cleaning up UI for slides without barcode
- Defining explicit location while adding new folders

## Bug fixes :

- Disabling of upload/download buttons on disconnection

From:

<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:

[https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.1.0](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.1.0)

Last update: **2021/06/18 17:25**

