

New features :

- Interfacing with the new My Pathomation API
- Adding create/rename/delete folder features for My Pathomation accounts
- Adding rename/delete slide features for My Pathomation accounts
- Adding drag & drop functionality for download
- Release of a self-containing package for Windows (No installation of JRE required)

General improvements :

- Detection of PMA.core version on connection
- Harmonization of upload/download UI/UX
- Cleaning up UI for slides without barcode
- Defining explicit location while adding new folders

Bug fixes :

- Disabling of upload/download buttons on disconnection

From:  
<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:  
[https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.1.0&rev=1590127613](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.1.0&rev=1590127613)

Last update: **2020/05/22 09:06**

