

## What's new in version 2.2.3?

PMA.transfer 2.2.3.9d6e296f92 is an improved version of the previous PMA.transfer with fixed bugs and new features in functionality.

### General improvements

- Upload slides and folders to the Amazon s3 cloud.
- Upload slides and folders recursively.
- Displaying information about slides in the cloud directories window.
- Adding a new folder in the cloud directory.

### New features in functionality

- Added the possibility to fill folders with slides and subfolders recursively, this means that you can load a folder containing subfolders containing both folders and slides.
- Added the possibility to unload multiple selected slides.
- Added the possibility to delete several selected slides from s3 the cloud directory at once.
- Added the possibility to create a new folder in the cloud S3 directory.

### Bug fixed :

- Bug fixed preventing more than one slide from being uploaded.

From:  
<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:  
[https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.2.3&rev=1662629305](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.2.3&rev=1662629305)

Last update: **2022/09/08 12:28**

