

# What's new in version 2.2.4.d67f7abd23?

PMA.transfer 2.2.4.d67f7abd23 is an improved version of the previous PMA.transfer with fixed bugs and new features in functionality.

## General improvements

- In this build version upload and download slides and folders execute whru Maven dependency(2.0.0.117)
- The filling of progressive bars executes thru the callback of the HTTP interceptor.
- Upload and download methods allow executing transfer slides to different storages: clouds(Azure and Amazon) and remote servers (file storage)

## New features in functionality

- Added the possibility to rename slide on local PC
- Added the possibility to rename folder on local PC

## Bug fixed

From:  
<https://docs.pathomation.com/PMA.transfer/> - PMA.transfer

Permanent link:  
[https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.2.4&rev=1677502628](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.2.4&rev=1677502628)

Last update: **2023/02/27 15:57**

