

# What's new in version 2.2.4.d67f7abd23?

PMA.transfer 2.2.4.d67f7abd23 is an improved version of the previous PMA.transfer with fixed bugs and new features in functionality.

## General improvements

- In this build version upload and download slides and folders execute whru Maven dependency(2.0.0.117)
- The filling of progressive bars executes thru the callback of the HTTP interceptor.
- Upload and download methods allow executing transfer slides to different storages: clouds(Azure and Amazon) and remote servers (file storage)

## New features in functionality

- Added the possibility to rename slide on local PC
- Added the possibility to rename folder on local PC

## Bug fixed

- Fixed progressive bar filling calculation system for slides with the multipart format as .vsi and .mrxs

## Slide formats

**Slide formats with PMA.transfer has been tested:**

```
* .mrxs
* .svs
* .ndpi
* .vsi
* .jpg
* .png
* .zif
* .tiff
```

Last update: 2023/02/27 16:05 what\_s\_new\_in\_version\_2.2.4 [https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.2.4&rev=1677503100](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.2.4&rev=1677503100)

---

From: <https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link: [https://docs.pathomation.com/PMA.transfer/doku.php?id=what\\_s\\_new\\_in\\_version\\_2.2.4&rev=1677503100](https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.2.4&rev=1677503100)

Last update: **2023/02/27 16:05**

