

What's new in version 2.2.4.d67f7abd23?

PMA.transfer 2.2.4.d67f7abd23 is an improved version of the previous PMA.transfer with fixed bugs and new features in functionality.

General improvements

- In this build version upload and download slides and folders execute whru Maven dependency(2.0.0.117)
- The filling of progressive bars executes thru the callback of the HTTP interceptor.
- Upload and download methods allow executing transfer slides to different storages: clouds(Azure and Amazon) and remote servers (file storage)

New features in functionality

- Added the possibility to rename slide on local PC
- Added the possibility to rename folder on local PC

Bug fixed

- Fixed progressive bar filling calculation system for slides with the multipart format as .vsi and .mrxs

Slide formats

Slide formats with PMA.transfer has been tested:

```
* .mrxs
* .svs
* .ndpi
* .vsi
* .jpg
* .png
* .zif
* .tiff
```

Last
update: 2023/02/27 16:05 what_s_new_in_version_2.2.4 https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.2.4&rev=1677503115

From:
<https://docs.pathomation.com/PMA.transfer/> - **PMA.transfer**

Permanent link:
https://docs.pathomation.com/PMA.transfer/doku.php?id=what_s_new_in_version_2.2.4&rev=1677503115

Last update: **2023/02/27 16:05**

