

# Glossary

- WSI - Whole Slide Image
- Annotation - a graphical object or layer of object added to a WSI
- Collection - a grouping of WSI that share a common context or purpose
- Slidebox - a grouping of multiple Collections and/or slides with added metadata for complex projects
- Snapshot - a picture of a specific area in a WSI
- Tags - short text describers for identifying WSI
- Free account - an personal account without costs, with restricted features
- My Workspace - a personal account with a subscription plan, offering advanced features
- Organization -a group account for multiple users
- Education and Collaboration - a subscription plan for a group of users with access to courses and quizzes features
- Case Sharing and AI - a subscription plan for a group of users with access to advanced case sharing and AI features

From:

<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:

<https://docs.pathomation.com/my/doku.php?id=glossary>

Last update: **2026/03/23 14:00**

