

# How to make a Course

The My Course option allows you to build content that consists of interactive images and text that you can share with students or collaborators. To create a new course, start by (1) +New course, once the course is created you can add/edit the content by opening the edit pane (2), the end result is shown when you click (3), you can simply share a course via a link like images or collections.

From:  
<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:  
[https://docs.pathomation.com/my/doku.php?id=how\\_to\\_make\\_a\\_course&rev=1677138059](https://docs.pathomation.com/my/doku.php?id=how_to_make_a_course&rev=1677138059)

Last update: **2023/02/23 10:40**

