

My Pathomation comes with a set of tools to annotate and draw on your slides have been updated. The available tools are listed here: 1/ Straight line & measurement tool 2/ Arrow & text tool 3/circle tool 4/ Square tool 5/ Freehand line and measurement tool 6/ Closed freehand tool 6/ text tool An illustration of the results is shown here

The annotations you make on a slide will persist if you include this slide in a collection, slide box or course

From:

<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:

https://docs.pathomation.com/my/doku.php?id=how_to_use_annotations&rev=1675354971

Last update: **2023/02/02 19:22**

