


## How to annotate slides?

My Pathomation comes with a set of tools to annotate and draw on your slides have been updated. The available tools are listed here: 1/ Straight line & measurement tool 2/ Arrow & text tool 3/circle tool 4/ Square tool 5/ Freehand line and measurement tool 6/ Closed freehand tool 6/ text tool An illustration of the results is shown here



The annotations you make on a slide will persist if you include this slide in a collection, slide box or course

From:  
<https://docs.pathomation.com/my/> - My Pathomation

Permanent link:  
[https://docs.pathomation.com/my/doku.php?id=how\\_to\\_use\\_annotations&rev=1675355849](https://docs.pathomation.com/my/doku.php?id=how_to_use_annotations&rev=1675355849)

Last update: 2023/02/02 19:37

