


Member View

As a manager or editor of an org, sometimes you want to test how your Quizzes work, or want to complete a Quiz made for you by another user.


To do this you can enter 'Member View'

First click your initials in the top right and click 'Enter Member View'



Now experience your content as your users do! When you're ready to return to manager view, click the button in the top right.

Your quiz results will be saved with your email address.



From:

<https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link:

https://docs.pathomation.com/my/doku.php?id=member_view&rev=1725186818

Last update: **2024/09/01 13:33**

