



# Annotations in collections

You can also edit a slide's annotations by clicking the "Edit annotations" button (2). These annotations will only be visible in the collection.



Collection specific annotations can be downloaded for use in alternate software. Start in edit and then select download annotations.



From: <https://docs.pathomation.com/my/> - My Pathomation

Permanent link: [https://docs.pathomation.com/my/doku.php?id=other\\_features\\_of\\_collections&rev=1728380371](https://docs.pathomation.com/my/doku.php?id=other_features_of_collections&rev=1728380371)

Last update: 2024/10/08 12:39

