

My Pathomation Start Page

My Pathomation is a cloud platform offering services for digital pathology applications.




Features

With My Pathomation you can store and share images for all scanner file formats. It comes with advanced features like annotation tools, collection organisation, access management and student course building. Registered users can store whole slide images (WSIs) and share them in an easy and universal way.

Whole slide images can be organised in folders but also be part of Collections, which in turn can be part of Slideboxes.

The tag features provides an additional layer of data for identification.



Organisations

Organisations in My Pathomation offer many interesting new features including featuring your fully interactive slides in **courses** and **quizzes**, and adding your own ****personal branding****.


With these powerful tools you can create your own content which integrate your fully interactive slides with course material, your own logos / branding, and many different question types including multiple choice, essay open field, number input and more.

Organisations have multiple user profiles including managers, content editors and members, each with different access rights.

My Pathomation also comes with a course builder tool and easy setup and management of groups and access right to all content.

Storage Options and Pricing


My Pathomation offers free image storage up to 10GB. For storing larger volumes of images, a monthly, volume-based subscription fee is required starting at 30 EUR for 25GB storage. Users can increase and decrease their required storage volume as required per organisation or personal workspace. There is no upper limit for the storage volume.



We Support You(r file types)

The My Pathomation cloud platform is compatible with all routine whole slide image formats. If you find a filetype that's not included that you'd like to use, we'd love to hear from you so we can integrate it into our platform!

BRIGHTFIELD	FLUORESCENT
3DHistech MRXS	3DHistech MRXS
Aperio / Leica CWS	Aperio AFI
Aperio LEICA SCN	Aperio LEICA SCN
Aperio LEICA SVS	Aperio LEICA SVS
GE Omnyx / Inspirata RTS	Hamamatsu NDPIS
GE Omnyx JP2	Nikon ND2
Generic JPG	Nikon TIFF
Generic PNG	Leica LIF
Generic TIFF	Olympus OIR
Hamamatsu NDPI	Olympus VSI
Hamamatsu VMS	Open Microscopy OME-TIFF
Huron Technologies TIFF	Perkin Elmer QPTIFF
JPEG XR encoded slides	Zeiss CZI
Jpeg XR JXR	Zeiss LSM
KFBio KFB	Zeiss ZVI
Menarini DSight Raw INI	
Microsoft Deep Zoom .DZI	Z-STACKS
Motic MDS	3DHistech MRXS
Motic MDSX	Aperio LEICA SVS
Nikon ND2	Aperio LEICA SCN
Philips iSyntax	Hamamatsu NDPIS
Philips TIFF	Huron Technologies TIFF
Olympus VSI	Leica LIF
Olympus Webview	Nikon ND2
Open Microscopy OME-TIFF	Objective Imaging (Glissando)
Objective Imaging (Glissando)	Open Microscopy OME-TIFF
Omero ZARR	Olympus OIR
Sakura SVSLIDE	Olympus VSI
Smartzoom SZI	Sakura SVSLIDE
Unic Tech TMAP	Zeiss ZVI
Ventana Roche BIF	Zeiss LSM
(with and without TIFF)	Zeiss CZI
Zeiss ZVI	
Zeiss CZI	
Zoomify ZIF	



From: <https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link: <https://docs.pathomation.com/my/doku.php?id=start&rev=1706608253>

Last update: **2024/01/30 12:50**

