


Your Dashboard landing page in My Pathomation

When you log in to you My Pathomation account you first see a 'My Dashboard' page grouping the key information of your My workspace or the organization you are logged in to. The dashboard:

1. has the name of the account you are logged in to, eg 'My workspace' or any of the organisations you have access to
2. shows key metrics about your account
3. lists all latest activities
4. allows you to click the items and then links directly to the relevant sections or slides
5. you can return to the dashboard at any time by clicking the logo



From: <https://docs.pathomation.com/my/> - **My Pathomation**

Permanent link: https://docs.pathomation.com/my/doku.php?id=your_dashboard_landing_page_in_my_pathomation

Last update: **2026/03/23 13:23**

