


Accessing the system

In order to access PMA.control, you must log in using your credentials. Depending on the configuration, self registration may be available, whereas you can create an account for yourself. In any other case, the credentials will be provided to you by the managers of the platform.



The screenshot shows the PMA.control login page. At the top left is the PMA.control logo. Below it is a 'Login' section with a text input field for 'Username, email' and a 'Password' field. A 'Keep me logged in' checkbox is checked. There is a 'Forgot password? Click here to register' link and a 'Sign In' button. At the bottom, the Pathomation logo is displayed with version information: 'v2.0.0.619' and 'PMA.core 2.0.0.622 ready'. Below the logo is the text 'Log in dialog'.

PMA.control is a web application designed to operate efficiently using modern web browsers. It is fully functional whether accessed from a desktop workstation or a touch enabled tablet. Devices with smaller screen size, such as smartphones, are compatible but not recommended, due to the limited available screen real estate.

From:
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:
https://docs.pathomation.com/pathotrainer/doku.php?id=accessing_the_system&rev=1596718357

Last update: **2020/08/06 15:52**

