


## Accessing the system

In order to access PMA.control, you must log in using your credentials.



The screenshot shows the PMA.control login page. At the top, there is a logo with a gear icon and the text 'PMA.control'. Below this is a 'Login' section with two input fields: 'Username, email' and 'Password'. The 'Password' field has a lock icon on the right. Below the password field is a checkbox labeled 'Keep me logged in' which is checked. To the left of the 'Sign In' button, there are links for 'Forgot password?' and 'Click here to register'. The 'Sign In' button is a dark grey rectangle with white text. At the bottom of the page, there is a logo for 'PATHOMATION' with the version 'v2.0.0.619' and 'PMA.core 2.0.0.622 ready' in green text. Below that is a link 'Log in dialog'.

PMA.control is a web application designed to operate efficiently using modern web browsers. It is fully functional whether accessed from a desktop workstation or a touch enabled tablet. Devices with smaller screen size, such as smartphones, are compatible but not recommended, due to the limited available screen real estate.

From:  
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:  
[https://docs.pathomation.com/pathotrainer/doku.php?id=accessing\\_the\\_system&rev=1597756760](https://docs.pathomation.com/pathotrainer/doku.php?id=accessing_the_system&rev=1597756760)

Last update: **2020/08/18 16:19**

