


Accessing the system

In order to access PMA.control, you must log in using your credentials.



The screenshot shows the PMA.control login page. At the top, there is a header with the PMA.control logo. Below the header, there is a 'Login' section with two input fields: 'Username, email' and 'Password'. The 'Password' field has a toggle icon for visibility. Below the password field, there is a checkbox labeled 'Keep me logged in' which is checked. To the right of the checkbox, there is a 'Sign In' button. Below the login fields, there is a link for 'Forgot password?' and a link for 'Click here to register'. At the bottom of the page, there is a logo for 'PATHOMATION' with the version 'v2.0.0.619' and 'PMA.core 2.0.0.622 ready' in green text. Below the logo, there is a link for 'Log in dialog'.

PMA.control is a web application designed to operate efficiently using modern web browsers. It is fully functional whether accessed from a desktop workstation or a touch enabled tablet.

Devices with smaller screen size, such as smartphones, are compatible but not recommended, due to the limited available screen real estate.

From:
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:
https://docs.pathomation.com/pathotrainer/doku.php?id=accessing_the_system&rev=1635697399

Last update: **2021/10/31 19:23**

