


# Accessing the system

In order to access Pathotrainer, you must log in using your credentials.



PMA.control

Login

Username, email

Password

Keep me logged in

Forgot password?  
[Click here to register](#)

Sign In

PATHOMATION  
v2.0.0.619  
PMA.core 2.0.0.622 ready

Log in dialog

Pathotrainer is a web application designed to operate efficiently using modern web browsers. It is fully functional whether accessed from a desktop workstation or a touch enabled tablet.

Devices with smaller screen size, such as smartphones, are compatible but not recommended, due to the limited available screen real estate.

From:  
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:  
[https://docs.pathomation.com/pathotrainer/doku.php?id=accessing\\_the\\_system&rev=1695384892](https://docs.pathomation.com/pathotrainer/doku.php?id=accessing_the_system&rev=1695384892)

Last update: **2023/09/22 15:14**

