


Accessing the system

In order to access Pathotrainer, you must log in using your credentials.



PMA.control

Login

Username, email

Password

☒ Keep me logged in

[Forgot password?](#) [Click here to register](#) [Sign In](#)

PATHOMATION
v2.0.0.619
PMA.core 2.0.0.622 ready
Log in dialog

Pathotrainer is a web application designed to operate efficiently using modern web browsers. It is fully functional whether accessed from a desktop workstation or a touch enabled tablet.

Devices with smaller screen size, such as smartphones, are compatible but not recommended, due to the limited available screen real estate.

From:
<https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link:
https://docs.pathomation.com/pathotrainer/doku.php?id=accessing_the_system&rev=1695384892

Last update: **2023/09/22 15:14**

