


# Adding Additional Material and Links

When [editing a project](#) the Additional material and Links tabs will allow you to add.... Additional material (i.e. non-WSIs) including videos or powerpoint presentations, and links to other websites.

## Additional Material

Navigate to the 'Additional material' tab, and press add




Now you can define:

1. The title of the content
2. Toggle if it's for supervisors' eyes only
3. Define the type: Introductory, Supplemental or Learning content
4. Specify the file for upload
5. Further organize the content by Topics which you can Add using the Add button, or select from a dropdown menu of topics you've already defined
6. Describe the additional material
7. And save when you're done

## Links

Navigate to the 'Links' tab, and press add



Now you can define:

1. The title you want to give to the link
2. The URL / web address you want to link to
3. Toggle if it's for supervisors' eyes only
4. Save when you're done!

From: <https://docs.pathomation.com/pathotrainer/> - **Pathotrainer**

Permanent link: [https://docs.pathomation.com/pathotrainer/doku.php?id=additional\\_material\\_and\\_links&rev=1701355218](https://docs.pathomation.com/pathotrainer/doku.php?id=additional_material_and_links&rev=1701355218)

Last update: **2023/11/30 17:40**

